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The cyanide ion-selective electrode may be used to measure the cyanide ion concentration in factory drainage, rivers, and lakes, or used in plant management in the pharmaceutical or plating industries.

In order to ensure correct usage of the ion electrode, be sure to read the operating instructions carefully before usage.

1 Contents of the Package

Cyanide Ion-selective Electrode 1
 Water proof abrasive paper (#1200) 1
 Instruction Manual 1

2 Specifications

| | |
|---|--|
| Concentration range | 0.03mg/L~2,600mg/LCN ⁻ (10 ⁻⁶ ~10 ⁻¹ mol/L CN ⁻) |
| pH range* ¹ | pH12~13 (2.6mg/LCN ⁻) |
| Temperature range | 0~50°C |
| Slope(mV/activity decade)* ² | -59mV |
| Response time T90* ³ | Less than 10 seconds |
| Material used for the body | Epoxy resin |
| Supported reference electrodes | |
| 8001-06T | 2080A-06T or 2535A-06T |
| 8001-10C | 2060A-10T or 2565A-10T |
| Reference electrode | 3.33mol/LKCl(Cat.No.300) |
| Type of internal filling solution | (248.3 g/L KCl) |

(※1)

This is the permitted pH range for the sample liquid ion concentration value given inside the brackets. (The pH range varies slightly depending on the ion concentration).

(※2)

This is the change in the electrode potential when the ion concentration becomes 10 times denser (at 25 degrees C.)

(※3)

The response time when the concentration is changed from 10⁻⁴ mol/L to 10⁻² mol/L.

3 The Electrode


Fig.1

4 Electrode Characteristics

1. Detection characteristic

Figure 2 shows the relationship between the cyanide ion concentration and the electrode potential.


Fig.2

2. Influence of pH

In acid, cyanide ions combine with the hydrogen to become hydrogen cyanide, making it impossible to measure using the electrode. Be sure to maintain the pH at 12 or higher in the standard solution and the test sample. Figure 3 shows the relationship between the pH and the electrode potential output for the cyanide ion electrode.


Fig.3

Relationship between pH and HCN in hydrogen cyanide

| pH | 6 | 7 | 8 | 9 | 10 | 11 | 12 | 13 |
|---|------|------|------|------|------|------|------|------|
| $\frac{\text{HCN}}{\text{HCN}+\text{CN}^-}$ % | 99.9 | 99.3 | 93.3 | 58.1 | 12.2 | 1.37 | 0.14 | 0.14 |

The HCN dissociation constant $K_a = 7.2 \times 10^{-10}$ (25°C)

3. Influence of coexisting ions

The cyanide ion electrode is influenced by various other ions. The influence of the coexisting ions can be shown in the following equation. K is the selectivity constant and is used as a standard for quantitatively indicating the influence of the coexisting ions.

$$E = E_0 - \frac{2.303RT}{F} \log \left(a_{\text{CN}^-} + \frac{1}{K} a_x \right)$$

Where a_{CN^-} is the cyanide ion concentration, a_x is the concentration of interfering ions. The effects of the other ions becomes smaller as K becomes larger, or when the ion concentration in the test sample increases. The selectivity constant K for various ions is listed in Table-1.

Table-1. Selectivity constants

| Interfering Ion | Selectivity constant (Limit capacity) |
|-------------------|---------------------------------------|
| $S_2O_3^{2-}$ | 1 |
| I^- | 0.1 |
| S^{2-}, MnO_4^- | Not possible |

Note : "Not possible" means that these ions cause a chemical change in the ion sensitive membrane.

4. Influence of temperature

The calibration curve of the electrode varies with temperature. You must maintain the standard solution and the sample solution at the same temperature during measurement.

5 Preparations

1. Preparation of the standard solution

Accurately weigh out 0.25 g of commercially available guaranteed reagent grade potassium cyanide and dissolve it in 0.1 mol/L sodium hydroxide (4 g/L NaOH). This solution is your 100 mg/L CN^- standard solution. Dilute 100 mL, 10 mL, and 1 mL of this solution with 1L of the 0.1 mol/L sodium hydroxide to produce 10 mg/L, 1 mg/L and 0.1 mg/L CN^- standard solutions.

2. Electrode conditioning

Prepare the electrode before taking measurements by soaking it for approximately 10 minutes in the 1 mg/L CN^- standard solution.

6 Measuring

1. Equipment

Assemble the electrode and ion meter (potentiometer) as shown in figure 4.


Fig.4

- Stir the solution at a constant speed (approximately 300 to 500 rpm) using the magnetic stirrer.
- Because of heat generated from the stirrer, the temperature of the solution is slowly rising, causing a drift in the electrode potential. If this is the case, lay a sheet of adiabatic material such as styrofoam between the stirrer and the beaker in order to prevent temperature changes.
- The stirring may generate noise in the electric potential of the electrode when measuring CN^- concentrations of greater than 100 mg/L. If this happens, and you are using an ion meter with a smoothing function, enable the smoothing function. An alternative is to dilute the test sample to less than 10 mg/L using 0.1 mol/L sodium hydroxide. This will reduce the effects of the stirring and allow stable measurements to be taken.
- Use a recorder if required.

2. Preparation of the sample solution

Add commercially available guaranteed reagent grade NaOH to the sample solution until it becomes 0.1 mol/L (4 g/L NaOH).

3. Operation

Create a calibration curve for the electrode using 2 or 3 types of standard solutions with 10-fold differences in concentration. This step is both for preparing the calibration curve, and also for verifying operation. The next step is to make quantitative measurements of the sample solutions using one of the following methods.

- Direct potentiometric Method
- Known Addition Method
- Known Subtraction Method
- Potentiometric Titration Method

Consult the JIS K 0122 standards (General Rules for Ion Electrode Methods) which covers each of these measurement techniques in detail. Here is a brief overview of (a) Direct Potentiometric Method.

With the Direct Potentiometric Method, you first create a calibration curve using two types of standard solutions with ion concentrations slightly higher and slightly lower than the ion concentration in the sample solution. You then measure the electrode potential in the sample solution, and arrive at the quantification value using the calibration curve.

7 Precautions

1. Do not rub or hit the ion sensitive membrane of the electrode. It may damage the membrane and degrade its special characteristics.
2. If direct sunlight falls on the ion sensitive membrane, the light may cause changes in the electrode potential. Use a brown beaker when using the equipment under such conditions.
3. Before soaking the electrode in a liquid with a different concentration, be sure to wash it first with deionized water.
4. The electrode connector requires a high degree of insulation. Take care not to touch this part with dirty hands or dampen it.

8 Troubleshooting

1. Check whether the problem is in the ion meter (potentiometer).
2. If the response from the electrode is slow, or its sensitivity has degraded, take the following steps to revive the ion sensitive membrane.
 - Drop a few drops of water onto the water proof abrasive paper included in the package. Lightly polish the response membrane by holding the electrode vertically over the paper and rubbing it in circles, then recalibrate the electrode.

9 Storage

- When storing the electrode for long periods, wash the ion sensitive membrane on the electrode with deionized water, and wipe it with some gauze. Fit the rubber cap included with the package over the response membrane and store the electrode in a cool dark place with a dry atmosphere.

HORIBA

Main Office/Factory: 2-Banchi, Kisshouinmiya-no-higashi-machi Minami-ku, Kyoto 601-8510
Tel: (075) 313-8121
Tokyo Sales Office: 1-7-8 Higashi-kanda Chiyoda-ku, Tokyo 101-0031
Tel: (03) 3861-8231